

# Detector Development for Nuclear Medical Imaging of Breast Cancer

**William W. Moses**  
Lawrence Berkeley National Laboratory  
January 12, 2000

## Outline:

- Breast Cancer / Mammography
- Positron Emission Tomography
- Detector Development at LBNL


## Breast Cancer Statistics

- 1 Woman in 8 Develops Breast Cancer
- 186,000 New Cases per Year in U.S.
- 46,000 Deaths per Year in U.S.
- Second Leading Cause of Cancer Death in Women
  - (slightly lower than lung cancer)

**BUT**

- Early Detection ⇒ Very High Cure Rate
  - Cancer hasn't spread into the axillary nodes ⇒ less likely to die of breast cancer than people who have never had breast cancer

## Early Detection via X-ray Mammograms



- Breast compressed to ~5 cm thickness.
- X-rays pass through breast, then “Bucky Grid” for scatter rejection.
- X-rays that pass through breast and Bucky Grid are detected in 2-D detector.
- Exposure controlled by integrated current in solid state (Si) detector.

## Low Energy X-Ray Image

## Mammography Images



**Abnormal Mass**




**Microcalcifications**



**Enlargement**

- Typical Contrast ~2%
- 25% – 80% False Positive Readings

## Advantage of Nuclear Medical Imaging



- X-Ray Mammography images **structure** (i.e. electron density)
- Nuclear Medicine images **function** (i.e. metabolic rate)

## Images the Tumor, Not Secondary Effects

## Positron Emitting Drug for Mammography:


### FDG (<sup>18</sup>F-fluoro-deoxy-glucose):

- >4x greater uptake in tumors than normal tissue
- >95\*% Sensitivity (i.e. probability that it has high uptake in breast cancer tumors)
- >80\*% Selectivity (i.e. probability that objects with high uptake are actually breast cancers).

\* for tumors  $\geq 1$  cm diameter

- Excellent Tracer for Imaging Breast Cancer
- Emits Back-to-Back 511 keV Photons


## PET Images of Breast Cancer Patient



- Tumors Easily Seen (~5 mm spatial resolution)
- Typical Contrast ~400%

## Positron Emission Tomography (PET)


Ring of Photon Detectors



- Patient injected with drug having  $\beta^+$  emitting isotope.
- Drug localizes in patient.
- Isotope decays, emitting  $\beta^+$ .
- $\beta^+$  annihilates with  $e^-$  from tissue, forming back-to-back 511 keV photon pair.
- 511 keV photon pairs detected via time coincidence.
- Positron lies on line defined by detector pair.

Forms Planar Image of a "Slice" Through Patient


## Multi-Layer PET Cameras



- Can image several slices simultaneously
- Can image cross-plane slices
- Can remove septa to increase efficiency ("3-D PET")

## Planar Images "Stacked" to Form 3-D Image


## Form Images Using Computed Tomography




**By measuring all 1-dimensional projections of a 2-dimensional object, you can reconstruct the object**

## Advantage of Computed Tomography

Planar X-Ray



Computed Tomography



Allows Separation of Objects on Different Planes



Images courtesy of Robert McGee, Ford Motor Company

## PET Cameras




General Electric




Siemens / CTI

## Conventional PET Detector Module




**Good Performance, Inexpensive, Easy to Pack**

## Crystal Identification with Anger Logic



**Can Decode Up To 64 Crystals with BGO**

## Positron Emission Mammography (PEM)




### Compared to 2D PET

- 30x Higher Sensitivity
- 2x Finer Resolution
- 4x Lower Dead Time
- 10x Lower Cost

**PET Camera Optimized for Breast Imaging**

## PET Camera for Axillary Nodes (ANPET)




### Compared to 2D PET

- 10x Higher Sensitivity
- 2x Finer Resolution
- 4x Lower Dead Time
- 10x Lower Cost

**Determine Whether Disease has Spread**




## Need Depth of Interaction Measurement



Required DOI Resolution is ~5 mm

## LBNL Detector Module Design



- PMT Provides Timing Pulse and Energy Discrimination
  - PD Array Identifies Crystal of Interaction
  - PD / (PD+PMT) Measures Depth of Interaction


## Lutetium Orthosilicate (LSO) Scintillator



	BGO	LSO
Formula	$\text{Bi}_4\text{Ge}_4\text{O}_{12}$	$\text{Lu}_2\text{SiO}_5:\text{Ce}$
Atten. Len.	1.1 cm	1.2 cm
Photo Fract.	43%	34%
Photon/MeV	8,200	25,000
Decay Time	300 ns	40 ns
Emission $\lambda$	480 nm	415 nm
Radioactive?	No	Yes (300 dps/cc)


## Recently Developed Scintillator

## Conventional PIN Photodiode Array




- Long Leads to Amplifier  $\Rightarrow$  Extra Capacitance  $\Rightarrow$  Extra Noise
- Vias at Edge  $\Rightarrow$  Dead Area

## LBNL Back-Side Illuminated Photodiode




- Short Leads  $\Rightarrow$  Low Capacitance  $\Rightarrow$  Less Noise
  - No Vias at Edge  $\Rightarrow$  Less Dead Area
  - Very Low Dark Current  $\Rightarrow$  Less Shot Noise

## Photodiode Array Performance



**LBNL Photodiode has Superior Signal & Noise**


## Amplifier Block Diagram



### Amplifier Characteristics


- Rise Time (0.2 – 100  $\mu$ s)
- Dark Current Sink (0 – 1000 pA)
- Gain (individual channel control, 50 – 200 mV / 1000  $e^-$  input)
- Fall Time (0.2 – 100  $\mu$ s)
- Power / Channel (1 – 20 mW)

## Amplifier Performance




**Plan to Operate at 0.5 – 1.0  $\mu$ s Shaping Time**

## The Crystal Identification Problem




**Provide Digital Address of Input with Highest Voltage in <100 ns**

## The "Winner Take All" Solution




- Amplifiers Fight for Limited Supply Current
- Output is Digital Bit Identifying Winner

## WTA Propagation Time



Well Under 100 ns

## 64 Channel IC Prototype Completed



### Inputs

- 64 PD Signals
- Sample & Hold
- Rise, Fall, & Reset Control

### Outputs


- Voltage of Winner
- Address of Winner

### Input / Output

- Serial Control (Xicor Protocol)


Submitting New Version Soon...

## Prototype Detector Module



## 64 Channel Prototype Built & Tested

### Depth Estimator $PD / (PD + PMT)$ is Linear with Excitation Depth




- Depth Resolution is 7–13 mm fwhm
- Distribution Similar in All Channels

## Prototype Detector Performance

Quantity	Conventional Module	This Module	Effect
Timing Resolution	~4 ns fwhm	0.75 ns fwhm	Improved Random Rejection
Energy Resolution	20% fwhm	25% fwhm	Degraded Scatter Rejection
Crystal Size	4.0 x 4.4 mm	3.0 x 3.0 mm	Improved Position Resolution
Correct ID Fraction	53% (Center 1) 66% (Center 5)	65% (Center 1) 91% (Center 5)	Improved Position Resolution
D.O.I. Resolution	none	7–13 mm fwhm	Reduced Radial Elongation

- **Outperforms Conventional Detector Module**
- **Further Improvements Expected...**

## Single Photon Variant



### Good Single Gamma Emitting Drugs Exist

- Technetium-99m Sestamibi
- Emits 141 keV gamma


### Collimate w/ Parallel Holes in Pb Detect with CsI:TI Scintillator / Photodiode Array

- ~4x more light than LSO
- CsI:TI not dense or fast for PET

- **Uses Similar Components as PEM / ANPET Module**
- **Planar (not Tomographic) Images Formed**
- **Radiopharmaceutical Widely Available**



## Components of the Detector Module



**Same Custom IC as PET Module**

– Signal level is ~50% higher


**Same PD Array as PET Module**

– Optimized for 550 nm

**Collimators are “Off the Shelf”**

- Readout IC Mounted to Back of Photodiodes
- Easy to Scale Up to a Full Camera

## <sup>99m</sup>Tc Spectrum in LBNL Detector



- Average Energy Resolution =  $10.7 \pm 0.6\%$  fwhm
- Conventional Camera (NaI:TI / PMT) Typically 9% fwhm

## **Conclusions**

### **Nuclear Medicine Can Image Breast Cancer:**

Excellent Accuracy for  $\geq 1$  cm Tumors

### **Specialized Cameras Under Development:**

Higher Resolution & Efficiency, Lower Cost  
Image Breast and Axillary Nodes  
At LBNL and Other Institutions

### **Development at LBNL:**

Major Components Developed (Photodiode, Custom IC)  
Prototype Detector Module Shows Excellent Performance  
Lots of Work Yet To Do (Postdoc? – wwmoses@lbl.gov)

## **Acknowledgements**

### **People:**

Mark Andreaco\*, Jean François Beche, Eric Beuville, Tom Budinger, Steve Derenzo, Greg Gruber, Matt Ho, Steve Holland, Michelle Huesman, Ron Huesman, Jenny Huber, Atul Joshi, Issy Kipnis, Brad Krieger, Tony Ko, Emmanuele Mandelli, Oren Milgrome, Chuck Melcher\*, Marzio Pedrali-Noy, Ron Nutt\*, Ola Pettersson, Amy Shustz, Patrick Virador, Nadine Wang, Matthew Weng

\*CTI PET Systems, Inc.

### **Funding:**

National Institutes of Health  
Grants R01-CA67911, P01-HL25840, & R01-NS29655  
State of California Breast Cancer Research Program  
Grant 1RB-0068  
U.S. Department of Energy  
Office of Biological and Environmental Research